

Anti-SARS-CoV-2 Antibody

Catalog No.	Application
K79d9	SARS-CoV-2 kit

Description	Cultured in vitro, without contamination of mouse antibodies
Source	Mammalian cell
Species reactivity	React with S protein
Isotype	Monoclonal IgG
Form	Clear and colorless liquid without flocculent precipitate
Purity	Purity > 98% , purified by Protein A/G chromatography
Buffer	1×PBS, pH7.4
	Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage instructions	Store at 2 - 8°C (35.6 – 46.4°F) for up to 12 months or at <-20 °C (-4°F) for up to 24 months